


# BlueShaper PRO<sup>®</sup>

|    | MFD  | TAPER | Ø    |
|----|------|-------|------|
| Z1 | 0,90 | 10-2% | 0,14 |
| Z2 | 0,90 | 10-2% | 0,17 |
| Z3 | 0,88 | 5-6%  | 0,19 |
| Z4 | 0,91 | 5-6%  | 0,25 |

The image displays four different sizes of BlueShaper PRO drill bits, labeled Z1 through Z4. Each bit is shown with its corresponding MFD (Major Flute Diameter) and TAPER values. The bits are color-coded: Z1 is purple, Z2 is white, Z3 is yellow, and Z4 is red. The Ø (diameter) values are also provided for each bit. The bits are arranged horizontally, showing their full length and the tapered cutting edge.

# BlueShaper PRO<sup>®</sup>



**KFile**  
#10

Glide Path



**Z1** **Z2**  
RPM: 500  
Torque: 4

Shaping



**Z3** **Z4**  
RPM: 500  
Torque: 4

Finishing



*Irrigate, recapitulate with a 10 file and re-irrigate after removing each rotary file*